

EverCross™

.035 PTA Balloon

Take crossing to a  
*higher level.*

The next generation of .035" PTA Balloons



EVERCROSS™


# Introducing EverCross™

the next generation in PTA Balloon technology from ev3.

Specially designed to optimize control, lesion entry, and trackability, EverCross delivers both superior performance and an unmatched range of balloon sizes.

## Expand .035" PTA options for your patients

- Unmatched 200 mm and 150 mm balloon lengths reduce unnecessary inflations in long lesions
- Balloons as short as 15 mm for treating focal lesions
- RBP up to 20 atms for treatment of challenging lesions





## Stay on the right track


- Strong catheter shaft design provides added pushability
- Flexible balloon design keeps you on track
- CrossCoat™ Hydrophilic Coating improves delivery and lesion crossing

## Get your point across

- CrossTec™ Beveled Tip Technology is specially designed for reliable lesion crossing
- Smallest tip entry profile available
- Ultra low crossing profile



LESION ENTRY PROFILE (in)\*



\* ev3 data on file

# EverCross™ .035" PTA Balloon

Product Code 135 cm catheter length	Product Code 90 cm catheter length	Product Code 40 cm catheter length	Balloon Diameter (mm)	Balloon Length (mm)	Rated Burst Pressure (atm)	Recommended Introducer Sheath (F)
AB35W03020135	AB35W03020080		3.0	20	20	5
AB35W03030135	AB35W03030080		3.0	30	20	5
AB35W03040135	AB35W03040080		3.0	40	20	5
AB35W03060135	AB35W03060080		3.0	60	20	5
AB35W03080135	AB35W03080080		3.0	80	20	5
AB35W03100135	AB35W03100080		3.0	100	20	5
AB35W03120135	AB35W03120080		3.0	120	20	5
AB35W03150135	AB35W03150080		3.0	150	20	5
AB35W03200135	AB35W03200080		3.0	200	20	5
AB35W04015135	AB35W04015080		4.0	15	20	5
AB35W04020135	AB35W04020080		4.0	20	20	5
AB35W04030135	AB35W04030080		4.0	30	20	5
AB35W04040135	AB35W04040080		4.0	40	20	5
AB35W04060135	AB35W04060080		4.0	60	20	5
AB35W04080135	AB35W04080080		4.0	80	20	5
AB35W04100135	AB35W04100080		4.0	100	20	5
AB35W04120135	AB35W04120080		4.0	120	20	5
AB35W04150135	AB35W04150080		4.0	150	20	5
AB35W04200135	AB35W04200080		4.0	200	20	5
AB35W05015135	AB35W05015080		5.0	15	18	5
AB35W05020135	AB35W05020080	AB35W05020040	5.0	20	18	5
AB35W05030135	AB35W05030080	AB35W05030040	5.0	30	18	5
AB35W05040135	AB35W05040080	AB35W05040040	5.0	40	18	5
AB35W05060135	AB35W05060080	AB35W05060040	5.0	60	18	5
AB35W05080135	AB35W05080080	AB35W05080040	5.0	80	18	5
AB35W05100135	AB35W05100080		5.0	100	18	5
AB35W05120135	AB35W05120080	AB35W05120040	5.0	120	18	5
AB35W05150135	AB35W05150080		5.0	150	18	5
AB35W05200135	AB35W05200080		5.0	200	18	5
AB35W06015135	AB35W06015080		6.0	15	14	5
AB35W06020135	AB35W06020080	AB35W06020040	6.0	20	14	5
AB35W06030135	AB35W06030080		6.0	30	14	5
AB35W06040135	AB35W06040080	AB35W06040040	6.0	40	14	5
AB35W06060135	AB35W06060080		6.0	60	14	5
AB35W06080135	AB35W06080080	AB35W06080040	6.0	80	14	5
AB35W06100135	AB35W06100080		6.0	100	14	5
AB35W06120135	AB35W06120080	AB35W06120040	6.0	120	14	5
AB35W06150135	AB35W06150080		6.0	150	14	5
AB35W06200135	AB35W06200080		6.0	200	12	6
AB35W07015135	AB35W07015080		7.0	15	14	5
AB35W07020135	AB35W07020080	AB35W07020040	7.0	20	14	5
AB35W07030135	AB35W07030080		7.0	30	14	5
AB35W07040135	AB35W07040080	AB35W07040040	7.0	40	14	5
AB35W07060135	AB35W07060080	AB35W07060040	7.0	60	14	5
AB35W07080135	AB35W07080080		7.0	80	14	5
AB35W07100135	AB35W07100080		7.0	100	14	6
AB35W07120135	AB35W07120080		7.0	120	14	6
AB35W07150135	AB35W07150080		7.0	150	14	6
AB35W07200135	AB35W07200080		7.0	200	14	6
AB35W08020135	AB35W08020080	AB35W08020040	8.0	20	14	6
AB35W08030135	AB35W08030080		8.0	30	14	6
AB35W08040135	AB35W08040080	AB35W08040040	8.0	40	14	6
AB35W08060135	AB35W08060080	AB35W08060040	8.0	60	14	6
AB35W08080135	AB35W08080080		8.0	80	14	6
AB35W09020135	AB35W09020080		9.0	20	12	6
AB35W09030135	AB35W09030080		9.0	30	12	6
AB35W09040135	AB35W09040080		9.0	40	12	6
AB35W09060135	AB35W09060080		9.0	60	12	6
AB35W09080135	AB35W09080080		9.0	80	12	6
AB35W10020135	AB35W10020080		10.0	20	11	6
AB35W10030135	AB35W10030080		10.0	30	11	6
AB35W10040135	AB35W10040080		10.0	40	11	6
AB35W10060135	AB35W10060080		10.0	60	11	6
AB35W12020135	AB35W12020080		12.0	20	10	7
AB35W12040135	AB35W12040080		12.0	40	10	7
AB35W12060135	AB35W12060080		12.0	60	10	7

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.  
CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

## NEUROVASCULAR | PERIPHERAL VASCULAR

Access · Balloons · Carotid · Embolic Coils · Embolic Protection · Liquid Embolics · Plaque Excision · Procedural Support · Retrieval Devices · Stents

**ev3 Corporate**  
World Headquarters  
Peripheral Vascular  
3033 Campus Drive  
Plymouth, MN 55441  
USA

**ev3 Neurovascular**  
9775 Toledo Way  
Irvine, CA 92618  
USA  
PH +1 949 837 3700  
FX +1 949 837 2044

**ev3 Europe**  
International Headquarters  
106-108 rue La Boétie  
75008 Paris  
France  
PH +33 156 88 59 10  
FX +33 156 88 59 11

PH +1 763 398 7000  
FX +1 763 398 7001  
Cust Srvc +1 800 716 6700

www.ev3.net

